

## DATASHEET

# Simplify clinical workflows and secure device access for vital signs devices

Fast, secure access to Baxter's Welch Allyn Connex vital signs devices with single sign-on powered by Imprivata

Your need to balance the benefits of vital signs monitor connectivity to your EHR with the need for strong user authentication doesn't need to mean sacrificing workflow simplicity for security. Imprivata Enterprise Access Management with single sign-on (formerly Imprivata OneSign) integrates with Baxter's Welch Allyn Connex vital signs monitors to enable fast, secure badge-tap access at the bedside.



## Protect patient data and save valuable time by extending badge access that clinicians are already using



### Streamline clinical workflows

Securely log in to devices without manually entering your username and password so you can spend more time focused on patient care.\* With the Welch Allyn dual function barcode scanner and RFID reader, clinicians have the ease of using one device to sign on to the vital signs monitor and to scan their patient's wristband for a completely connected workflow.




### Control device access

A simple badge scan or tap provides a quick authentication at the bedside, granting device access only to authorized users.

## A connected clinical workflow with single sign-on

The process is simple. Users start by opening a session with Imprivata – either at an EHR-enabled workstation, on a vital signs monitor, or elsewhere. This authentication is then extended to a deployed Connex vital signs device. Access is granted in a single step, allowing clinicians to take patient vitals and send data to the EHR.



**Contact your Imprivata representative or visit [imprivata.com](https://imprivata.com) to learn more. Work with your Imprivata representative to validate usability with your specific RFID type.**

**Rx only.** For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

Note: RFID support on Connex Spot Monitor software versions 1.32 or higher and Connex Vital Signs Monitor 2.42 or higher.

\*You must first sign in with Imprivata in accordance with your facility's policies. Then you can simply tap or scan your badge on select Connex devices to log in to the device.



**Baxter.com | Baxter International Inc.**

Baxter, Connex and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.



Imprivata delivers simple and secure access management solutions for healthcare and other mission-critical industries to ensure every second of crucial work is both frictionless and secure. The Imprivata platform of innovative, interoperable access management and privileged access security solutions enables organizations to fully manage and secure all enterprise and third-party identities to facilitate seamless user access, protect against internal and external security threats, and reduce total cost of ownership.

For more information, please contact us at:

**Global headquarters USA**

Waltham, MA

**Phone:** +1 877 663 7446

[www.imprivata.com](http://www.imprivata.com)

**European headquarters**

Uxbridge, England

**Phone:** +44 (0) 208 744 6500

[www.imprivata.com/uk](http://www.imprivata.com/uk)

**Germany**

Langenfeld

**Phone:** +49 (0) 2173 99 385 0

[www.imprivata.com/de](http://www.imprivata.com/de)

**Australia**

Melbourne

**Phone:** +61 3 8844 5533

Copyright © 2025 Imprivata, Inc. All rights reserved. Imprivata is a registered trademark of Imprivata, Inc. in the U.S. and other countries. All other trademarks are the property of their respective owners.