

# Palm scanner care guide

## Taking care of your palm scanner

The palm scanner, while tolerant of residues, performs optimally when kept clean!

Carefully follow the recommendations below to ensure optimal use.

## Procedures for cleaning

### Cleaning instructions

1. Allow the patient to see you clean the scanner
2. Use alcohol wipes on any part of scanner
  - Note: **Do not** use CaviWipes™ (due to composition of the cleaner)
3. Use a microfiber towel to clean the sensor glass
4. **Do not** use liquids to clean the sensor glass
5. **Do not** take the scanner apart. Notify Support if it is not functioning properly
6. Clean the plastic scanner housing of the scanner whenever you clean the other non-critical health care equipment in your area
7. Use of a 1:10 bleach solution is acceptable on the plastic housing only, though it may accelerate wear on the Imprivata logo